DetectairTM
technologies

HERMES
Lab Systems

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NATIONAL INSTITUTE OF CHILDREN'S DISEASES
Limbová 1, 833 40 Bratislava, Slovak Republic

On April 3, 2023, we also took samples to determine the microbial settlement of the environment in operating theatres 2 and 3, immediately after measuring the air.

Naturally the following measures were adopted:

- regular disinfection of all surfaces and spaces using disinfectants according to the disinfection programme
- germicidal irradiation according to regulations.

There was also the apparent use of aerosol disinfection according to the key mentioned above.

SN	Assessed area	Primoculture	After multiplication	Vyhodnotenie
1.	OT 2 – rubber keyboard	POS	POS	Negative
2.	OT 2 – operating table	POS	POS	Negative
3.	OT 2 – injection machine	POS	POS	Negative
4.	OT 2 – the floor	POS	POS	Negative
5.	OT 2 – lateral vent, air conditioning	POS	POS	Negative
6.	OT 2 – disinfection nozzle	POS	POS	Negative
7.	OT 2 – aerosol disinfection pipeline	POS	POS	Negative
8.	OT 2 – suction machine	POS	POS	Negative
9.	OT 2 – anaesthesia table, monitor	POS	POS	Negative
10.	OT 2 – anaesthesia table, surface	POS	POS	Negative
11.	OT 3 – instrument table	POS	POS	Negative
12.	OT 3 – operating table	POS	POS	Negative
13.	OT 3 – lateral vent, air conditioning	POS	POS	Negative
14.	OT 3 – ceiling vent, air conditioning	POS	POS	Negative
15.	OT 3 – headlight, operating lamp	POS	POS	Negative
16.	OT 3 – operating table, controller	POS	POS	Negative
17.	OT 3 – inside door handle	POS	Staphylococcus sp. – coagulase negat.	Negative
18.	OT 3 – disinfection nozzle	POS	POS	Negative
19.	OT 3 – the floor	POS	POS	Negative
20.	OT 3 – anaesthesia table, coloured	POS	POS	Negative

Evaluation of microbial contamination of environment

The results of the evaluation of microbial contamination of the environment are in 19 out of 20 cases evaluated as: "the soils remained sterile." The door handle was the only place where a microorganism, which is commonly found on hands, was found, however I cannot rule out my hand touching the handle before the actual sampling.

Recommendation

Based on the above-mentioned results, I recommend using aerosol disinfection as an additional measure to increase the safety of patients when providing healthcare. In the course of summer 2023 I will continue in the use and further testing of aerosol disinfection in operating theatres.

The evaluation was done on May 18, 2023 by

Mgr. Michal Adamišín, PhD.
Public health worker – sanitarian

Tel.: +421 (0)2 59371405
ID number: 00607231

Fax: +421 (0)2 59371884
TIN: 2020848368

E-mail: michal.adamisin@nudzch.eu
VAT number: SK 2020848368

NATIONAL INSTITUTE OF CHILDREN'S DISEASES
Limbová 1, 833 40 Bratislava, Slovak Republic

Results of microbial contamination of air and environment

Assessed workplace: Main operating theatres
Present: Mgr. Michal Adamišín, PhD. – public health worker/sanitarian
Head nurse, department nurse
Sampling date: Multiple dates

Reason for sampling:
Sampling to detect microbial contamination of air and environment for the further testing of aerosol disinfection provided in operating theatres.

Sampling date	Assessed workplace	Aerobic cultivation	Evaluation
19.9.2022	Operating theatre 2 (100 litres of air, 1 minute)	<i>Micrococcus luteus</i> (90 colonies in 100 litres of air) <i>Staphylococcus</i> sp. – coagulase negat. (25 colonies in 100 litres of air)	Negative
21.9.2022	Operating theatre 3 (100 litres of air, 1 minute)	<i>Micrococcus luteus</i> (90 colonies in 100 litres of air) <i>Staphylococcus</i> sp. – coagulase negat. (80 colonies in 100 litres of air)	Negative
3.4.2023	Operating theatre 2 (100 litres of air, 1 minute)	POS	Negative
	Operating theatre 3 (100 litres of air, 1 minute)	Aerial saprophytes (18 colonies in 100 litres of air) <i>Staphylococcus</i> sp. – coagulase negat. (25 colonies in 100 litres of air)	Negative

Legend:
POS: soils remained sterile, satisfactory sampling
Orange tint: positive sampling, unsatisfactory sampling
Green tint: soils remained sterile
The risk microorganism is marker in **bold**.

Methodology of measurement

For air sampling we used sampling sets with Sabouraud and blood agar, which we let air blow through by sucking it with an aeroskop for 1 minute with a volume of 100 litres of air, due to the environment being operating theatres.

On September 19 and 21, 2022 we had not yet been using the aerosol disinfection. We chose operating theatres 2 and 3 due to the distance of all tubes from the control panel being almost identical.

On April 3, 2023, we carried out an aerosol disinfection of the environment for 120 to 180 seconds every working day for almost three months before proceeding with sampling. We stopped using aerosol disinfection in operating theatre 3, 2 weeks before sampling; previously it had been used every working day, also for 120 to 180 seconds.

Evaluation of air measurement results

Above all, I will focus on the result from April 3, 2023, obtained from operating theatre 2, when neither pathogenic nor non-pathogenic microorganisms were identified by measuring the air, and thus we identified the result as "the soils remained sterile."

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Date of issue: 30.01.2023
Dátum vydania: -

Tested and elaborated by: Patrícia Vavrinec
test technician

Checked and approved by: Ing. Mario Zímko
Head of machinery and production
devices Testing Body

The results introduced in this test report are related to the test subject only and do not substitute other documents required by state supervisory authorities and according to other specific regulations. Test report can be reproduced or published as a whole, in parts only, with written approval of TSU test body.

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TEST REPORT
no. 220500317/EMC

Test name: Electromagnetic compatibility tests
Product name: HYDROFOGG

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Test - procedure method: see chap. 1
Date of test performance: see chap. 4
Test equipment: see chap. 5
Distribution: 1x - manufacturer
1x - TSU Piešťany, s.p.
Test conditions: see data mention in particular articles

Date of issue: 26.01.2023

Ing. Lukáš Vando
Technical head of RED and EMC
testing body

Martin Ryšly
Test engineer

The results introduced in this test report are related to the test subject only. Test report can be reproduced or published as a whole, in parts only, with written approval of TSU test body.

ELBACERT
CERTIFICATION BODY

IAF
INTERNATIONAL ACCREDITED

SNAS
Reg. No. 1939/006

CERTIFICATE
Certification body ELBACERT, joint-stock company
certifies by this that company

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has implemented and maintained
an **Environmental Management System**
in accordance with requirements of standard

ISO 14001:2015
for the scope:
Design, manufacture, sale, service of equipment for disinfection
and
dedusting of premises.

Certificate issue date: June 05, 2024
Certificate validity within: August 15, 2026
Certificate no.: 2024272

Mr. Alois Švajda
Certification body director

The certificate is granted on the basis of ISO Decision No. 23/2024.
The relevance of this certificate depends on the validity of certificate 2024272 in its validity.
The certificate remains valid subject to the periodic nature of the annual surveillance audits.

ELBACERT, joint-stock company, Kálna 5, 020 01 Žilina, Slovak republic
www.elbcert.com

BeeIQ

ELBACERT
CERTIFICATION BODY

IAF
INTERNATIONAL ACCREDITED

SNAS
Reg. No. 1939/006

CERTIFICATE
Certification body ELBACERT, joint-stock company
certifies by this that company

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has implemented and maintained
an **Quality Management System**
in accordance with requirements of standard

ISO 9001:2015
for the scope:
Design, manufacture, sale, service of equipment for disinfection
and
dedusting of premises.

Certificate issue date: June 05, 2024
Certificate validity within: August 15, 2026
Certificate no.: 2024272

Mr. Alois Švajda
Certification body director

The certificate is granted on the basis of ISO Decision No. 23/2024.
The relevance of this certificate depends on the validity of certificate 2024272 in its validity.
The certificate remains valid subject to the periodic nature of the annual surveillance audits.

ELBACERT, joint-stock company, Kálna 5, 020 01 Žilina, Slovak republic
www.elbcert.com

BeeIQ

ELBACERT
CERTIFICATION BODY

IAF
INTERNATIONAL ACCREDITED

SNAS
Reg. No. 1939/006

CERTIFICATE
Certification body ELBACERT, joint-stock company
certifies by this that company

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has implemented and maintained
an **Occupational Health and Safety Management System**
in accordance with requirements of standard

ISO 45001:2018
for the scope:
Design, manufacture, sale, service of equipment for disinfection
and
dedusting of premises.

Certificate issue date: June 05, 2024
Certificate validity within: August 15, 2026
Certificate no.: 2024274

Mr. Alois Švajda
Certification body director

The certificate is granted on the basis of ISO Decision No. 23/2024.
The relevance of this certificate depends on the validity of certificate 2024274 in its validity.
The certificate remains valid subject to the periodic nature of the annual surveillance audits.

ELBACERT, joint-stock company, Kálna 5, 020 01 Žilina, Slovak republic
www.elbcert.com

BeeIQ

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2025 © Všetky práva vyhradené.